

Camera Solutions



Camera Solutions for Olympus Trinocular Microscopes (1)







Upright Microscope

Applicable Model:

CX23T, CX33, CX43, BX43, BX46, BX53, BX63, SZX16, SZX10, SZX7, SZ61TR



	Solution 1	Solution 2	Solution 3
Deeply customized Android OS	✓	✓	
Pre-installed MS Office	✓	✓	
15.6-inch High Color Gamut Display	✓		
Output Methods			
5G WiFi	✓	✓	✓
USB			✓
HDMI	✓	✓	✓
WAN (Ethernet)			✓

Camera Solutions for Olympus Trinocular Microscopes (2)



Research Upright Microscope

Applicable Model:

BX43, BX46, BX53, BX63, etc.



	Solution 4	Solution 5
Deeply customized Android OS		
Pre-installed MS Office		
15.6-inch High Color Gamut Display		
Output Methods		
5G WiFi	✓	\checkmark
USB		\checkmark
HDMI		
WAN (Ethernet)	✓	✓


Camera Solutions for Olympus Trinocular Microscopes (3)



Inverted Microscope

Applicable Model:

CKX53, GX53, etc.



Functions Available:

	Solution 1	Solution 2	Solution 3
Deeply customized Android OS	✓	✓	
Pre-installed MS Office	✓	✓	
15.6-inch High Color Gamut Display	✓		
Output Methods			
5G WiFi	✓	✓	✓
USB			✓
HDMI	✓	✓	✓
WAN (Ethernet)			✓

① HDMI output is unavailable for 20MP cameras.

Camera Solutions for Nikon Trinocular Microscopes (1)



Upright Microscope

Applicable Model:

ECLIPSE Ei, Si, Ni, Ci, LV100, LV150, L200, L300, SMZ series



	Solution 1	Solution 2	Solution 3
Deeply customized Android OS	✓	✓	
Pre-installed MS Office	✓	✓	
15.6-inch High Color Gamut Display	✓		
Output Methods			
5G WiFi	✓	✓	\checkmark
USB			\checkmark
HDMI	✓	✓	\checkmark
WAN (Ethernet)			\checkmark

Camera Solutions for Nikon Trinocular Microscopes (2)



Research Upright Microscope

Applicable Model:

ECLIPSE Ni, Ci, etc.

Solution 4
Multi-viewing Workstation



Dovetail mount with 0.43X relay lens or C mount



Dual-sensor Camera



Dovetail mount with 0.63X lens

Functions Available:

	Solution 4	Solution 5
Deeply customized Android OS		
Pre-installed MS Office		
15.6-inch High Color Gamut Display		
Output Methods		
5G WiFi	✓	✓
USB		✓
HDMI		
WAN (Ethernet)	✓	✓

• HDMI output is unavailable for 20MP cameras.

Camera Solutions for Nikon Trinocular Microscopes (3)



Inverted Microscope

Applicable Model:

ECLIPSE Ts2, MA100N



Dedicated Multi-output Camera



Dovetail mount with 0.43X relay lens or C mount







with 0.43X relay lens



	Solution 6	Solution 7
Deeply customized Android OS		
Pre-installed MS Office		
15.6-inch High Color Gamut Display		
Output Methods		
5G WiFi	✓	✓
USB	✓	\checkmark
HDMI	✓	\checkmark
WAN (Ethernet)	✓	\checkmark

Camera Solutions for Leica Trinocular Microscopes (1)





Upright Microscope

Applicable Model:

DM500, DM750, DM2500, DM3500, DM3500, DM2700P, DM2700M, DM750P, etc.










	Solution 1	Solution 2	Solution 3
Deeply customized Android OS	✓	✓	
Pre-installed MS Office	✓	✓	
15.6-inch High Color Gamut Display	✓		
Output Methods			
5G WiFi	\checkmark	✓	✓
USB			✓
HDMI	✓	✓	✓
WAN (Ethernet)			✓

Camera Solutions for Leica Trinocular Microscopes (2)



Functions Available:


Dovetail mount with 0.43X relay lens or C mount

	Solution 4	Solution 5
Deeply customized Android OS		
Pre-installed MS Office		
15.6-inch High Color Gamut Display		
Output Methods		
5G WiFi	✓	✓
USB		✓
HDMI		
WAN (Ethernet)	✓	✓

• HDMI output is unavailable for 20MP cameras.

Dovetail mount with 0.63X lens

Camera Solutions for Leica Trinocular Microscopes (3)




Functions Available:

	Solution 2	Solution 3
Deeply customized Android OS	✓	
Pre-installed MS Office	✓	
15.6-inch High Color Gamut Display		
Output Methods		
5G WiFi	✓	✓
USB		✓
HDMI	✓	✓
WAN (Ethernet)		✓

• HDMI output is unavailable for 20MP cameras.

Camera Solutions for Leica Trinocular Microscopes (4)



with dedicated adapter

Stereo Microscope

Dedicated for stereo microscope

Functions Available:

Solution 6

All-in-one Smart Camera

with dedicated adapter

	Solution 6	Solution 7	Solution 8	Solution 9
Deeply customized Android OS	✓	✓		✓
Pre-installed MS Office	✓	✓		✓
15.6-inch High Color Gamut Display	✓			
Output Methods				
5G WiFi	✓	✓	✓	\checkmark
USB			✓	
HDMI	✓	✓	✓	✓
WAN (Ethernet)			✓	

with dedicated adapter

HDMI output is unavailable for 20MP cameras.

Camera Solutions for Zeiss Trinocular Microscopes (1)







Upright Microscope

Applicable Model:

Axio Examiner, Primo, Axio Imager, Axiolab, Axioscope, Axio Observer, Stemi 508, Axiolab 5, etc.

Solution 3







	Solution 1	Solution 2	Solution 3
Deeply customized Android OS	\checkmark	✓	
Pre-installed MS Office	\checkmark	✓	
15.6-inch High Color Gamut Display	✓		
Output Methods			
5G WiFi	\checkmark	✓	\checkmark
USB			\checkmark
HDMI	\checkmark	✓	\checkmark
WAN (Ethernet)			\checkmark

Camera Solutions for Zeiss Trinocular Microscopes (2)



Research Upright Microscope

Applicable Model:


Axio Examiner, Primo, Axio Imager, Axiolab, Axioscope, etc.

Solution 4
Multi-viewing Workstation



Dovetail mount with 0.43X relay lens or C mount







Dovetail mount with 0.63X lens

	Solution 4	Solution 5
Deeply customized Android OS		
Pre-installed MS Office		
15.6-inch High Color Gamut Display		
Output Methods		
5G WiFi	\checkmark	✓
USB		✓
HDMI		
WAN (Ethernet)	\checkmark	✓

[•] HDMI output is unavailable for 20MP cameras.

Camera Solutions for Zeiss Trinocular Microscopes (3)



Multi-viewing Workstation

Solution 2



Dovetail mount with 0.43X relay lens or C mount

Solution 3

Dual-sensor Camera



Dovetail mount with 0.43X lens

Functions Available:

	Solution 2	Solution 3
Deeply customized Android OS	✓	
Pre-installed MS Office	✓	
15.6-inch High Color Gamut Display		
Output Methods		
5G WiFi	✓	✓
USB		✓
HDMI	\checkmark	✓
WAN (Ethernet)		✓

HDMI output is unavailable for 20MP cameras.



LANOPTIK TECHNOLOGIES LTD

No. 72 Hongjing Street, Baiyun, Guangzhou, China. 510400
Phone: +86 20 3898 6017 | Fax: +86 20 3847 6076

LANOPTIK
Technologies

Website: http://www.lanoptik.com | Email: info@lanoptik.com